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ANALYSIS TECHNOLOGY

Malvern

Applications Report on Rheological Measurements of

CMC + CARBOPOL + PAPAINA
CMC + CARBOPOL + PAPAINA (different preparation)

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
Annex 2

Procedure

Exactly the same methods, equipment, geometry, etc. used for the original samples were used for this study - It is worth mentioning that it can be seen on plain view a brown coloring of the sample compared to the three previous test samples, and our lab assumes that it is due to rusting thereof.

Findings

Flow Curve


The yield stress of the original sample (blue line), is higher than the final sample sent by the client. With this test we want to show that possibly it has a higher stability, because at this low speeds (shearing) a higher viscosity of the product is needed in order not to have the separation or sedimentation effects, etc.

It is worth mentioning that under the application area (shear rates of approx $10^1 s^{-1}$), the viscosity of both samples is equal, which is evidence of a similar spreading trait.


Frequency studies

The two samples show a lineal viscoelastic region (LVER) high enough to show a good trait of the gel.

Frequency

There can be seen similar behaviors in its angle features of the phase angle in the frequency oscillatory sweep, but the original sample has a lower angle towards the frequencies that simulate the application (meaning: The angle of the phase is moved between 0° to 90°, being 0° something completely elastic, or solid, and 90° something completely viscous. For an application a higher viscosity is better at the point of "structure breakage")

Applications Report on Rheological Measurements of

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HUMAX PHARMACEUTICAL S.A.